

Investigation of Suspected and Unknown Micropollutants and Transformation Products from a Waste Water Treatment Plant with Full Scale Ozonation

Christoph Portner,¹ Olaf Scheibner,² Sebastian Westrup,³ Jochen Tuerk¹

¹Institut für Energie- und Umweltechnik e. V., IUTA (Institute of Energy and Environmental Technology), Duisburg, Germany; ²Thermo Fisher Scientific, Bremen, Germany, ³ Thermo Fisher Scientific, Dreieich, Germany



Overview

- Investigation of micropollutants and transformation products
- Waste water treatment with ozonation as final step for elimination
- Comparison of influent, effluent after conventional biological cleaning and effluent after full scale ozonation
- Combination of suspected and non target screening with LC-HRMS
- Identification of transformation products using *in silico* transformation and fragmentation

Introduction

Due to the extensive use of industrial chemicals, pesticides, personal care products and pharmaceuticals, many anthropogenic organic micropollutants and their biological metabolites have become ubiquitously detectable in the aquatic environment. The emission of these substances into the aquatic environment occurs through both direct and indirect discharge into the sewage in combination with insufficient biological or chemical degradation efficiency of conventional municipal waste water treatment plants (WWTPs). Against this background, samples from the influent, the effluent after the biological treatment and after the final ozonation from the WWTP in Duisburg-Vierlinden, Germany were picked up and analyzed using liquid chromatography high-resolution mass spectrometry (LC-HRMS). The data investigation was focused on occurrence of suspected target and non target micropollutants, biological metabolites and transformation resulting from the chemical treatment step.

Methods

Sample Preparation

Samples were taken from the WWTP influent, the effluent after the biological treatment and effluent after the ozonation. The samples were filtered and extracted by solid phase extraction (Oasis® HLB, Waters®).

Liquid Chromatography

For chromatographic separation an ultra high pressure system (Thermo Scientific Aria Transcend™) was used. A sample volume of 60 µL was injected onto a 100x2.1 mm analytic column (Thermo Scientific Hypersil Gold™ aQ). A 7 minute solvent gradient (Eluent A: Water + 0.1% formic acid, B: Methanol + 0.1% formic acid) from 5-99% B was applied resulting in a total cycle time of 15 minutes for chromatographic separation at a flow rate of 400µl/min. The column oven were set to 25 °C.

Mass Spectrometry


For mass spectrometric detection a high resolution mass spectrometry system (Thermo Scientific Exactive Plus™) with positive and negative electrospray ionization was used and run in permanently alternated mode of full scan and all ion fragmentation (AIF). A resolution setting of 70.000 (FWHM @ m/z 200) was used. A mass range of m/z 100 to 1500 was applied (resp. m/z 50 to 750 and resolution setting 70.000 FWHM for the AIF scans) to be prepared for all possible contaminants. The mass axis of the system was calibrated with the standard calibration mix once prior measurement.

Data Analysis

The data were analyzed in a widely automated workflow using Thermo Scientific TraceFinder™ 3.1, Thermo Scientific SIEVE™ 2.1 and Thermo Scientific Compound Discoverer™ 1.0 software.

Results

FISH scoring for Metoprolol with fragmentation




The FISH scoring takes into account all fragment information available in the raw data and combines it with *in silico* fragmentation data of the compounds used for screening. The figure shows such an annotated product ion spectrum of Metoprolol. Besides isotopic pattern matching for all found hits the FISH scoring serves for an additional confirmation of the obtained results.

Non target screening


For unknown screening, elimination of background signals and un-significant peaks is a major issue. To achieve this, a multiple mass defect filtering step runs before the component extraction. The multiple mass defect filters are connected to the components used for suspect screening, so a large quantity of signals is eliminated which is not likely to be related to the parent components of the study. The resulting masses can be used for further identification using online databases like ChemSpider™ or mzCloud™.

Comparison of transformation products of carbamazepine



The transformation products of carbamazepine in the influent (A), the effluent (B) and after the ozonation indicates the complex reaction pathways during the WWTP process.

Identification of transformation products during ozonation



Volcano plot of the samples effluent after the biological treatment vs. after the ozonation. The figure shows the compounds formed during the ozonation and the eliminated ones. Many peaks were found in both samples.

Database file name	Formula	<i>m/z</i> (Measured)	Mass error (ppm)	Fit isotopic pattern	Fit RT	Identified <i>in silico</i>
TP 379	C ₁₅ H ₁₀ N ₂ O ₂	251.0814	-0.3328	Yes	Yes	Yes
TP 377	C ₁₅ H ₁₂ N ₂ O ₃	269.0922	0.1911	Yes	n.a.	Yes
TP 386	C ₁₅ H ₁₄ N ₂ O ₂	255.1127	-0.2695	Yes	Yes	Yes

Identified transformation products of carbamazepine in the effluent of the full scale ozonation for the suspected target screening using a custom made database with 354 transformation products of 32 common micro-pollutants based on a literature research and own experiments. All compounds could be verified by the *in silico* screening approach.

Conclusions

- ▶ **Successful elimination of micropollutants with ozonation as final waste water treatment step**
- ▶ **Identification of transformation products using a combination of suspected target and non target screening**
- ▶ **New screening approach for the investigation of transformation pathways during the biological waste water treatment process and ozonation by *in silico* transformation and fragmentation**

Acknowledgement

The authors would like to thank Thorsten Bernsmann (CVUA-MEL) and Claudia Martins (Thermo Fisher Scientific) for their support and the opportunity to use the LC-HRMS system at CVUA-MEL.



www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. Sigma Aldrich is a registered trademark of Sigma-Aldrich Co. LLC. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific,
San Jose, CA USA is
ISO 9001:2008 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00
Russia/CIS +43 1 333 50 34 0

Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand