

A Novel UHPLC Orbitrap HRMS Approach for Monitoring Sulfur Fumigation Abuse in Drying Process of Traditional Chinese Medicine

Zhe Zhou,¹ Min Yang²

¹Thermo Fisher Scientific, Shanghai, China; ²Shanghai Institute of Materia Medica, Shanghai, China

Overview

Purpose: To demonstrate the application of fine isotope pattern filter with ultra high resolution mass spectrometry data for sulfur fumigation abuse monitoring in Traditional Chinese Medicine.

Methods: Sulfur fumigation processed herb powder was analyzed by Thermo Scientific™ Orbitrap Velos Pro™ high resolution mass spectrometer. MS data processing was performed on Thermo Scientific™ Compound Discoverer™ software.

Results: The fine isotope pattern screening method was successfully applied in study of sulfur-fumigated *Pueraria* species with twelve major sulfur derivatives, eleven of which were tentatively identified according to their HRMSⁿ data which is not known to occur naturally in these samples. The factors which may effect the accuracy of screening result were also evaluated.

Introduction

Traditional Chinese medicine has been used for disease treatment for thousands of years and still plays a very important role in current medical system in China and worldwide. Recently, due to the benefits of reducing the drying duration, preventing against insects, maintaining better appearance and even keeping moist to gain more weight, sulfur fumigation has been commonly employed for fast drying medicinal materials during post-harvest handling instead of traditional natural drying process. However, determination of SO₂ involves complicated reactions to release and capture free SO₂ and was not sensitive enough. Benefit from the ultra high resolving power of Orbitrap Velos Pro mass spectrometer up to 100,000 FWHM, a novel fine isotopic pattern screening method was developed according to the separation of ³⁴S with ¹³C2 and ¹⁸O isotopic pattern with mass shift of only 0.00981 amu.

Methods

Sample Preparation

Powdered samples (40 mesh, 0.2 g of Yege (YG, root of *Pueraria lobata* (Willd.) Ohwi), 0.1 g of Gegenqinlian tablets and Yufengningxin tablets, 0.8 g of Fenge (FG, root of *Pueraria thomsonii* Benth.) were extracted with 30 mL 30% ethanol for 30 minutes of ultrasonic extraction, filtered.

Liquid Chromatography:

Chromatographic separation was performed using an Thermo Scientific™ UltiMate 3000 Binary RSLC system, equipped with a binary solvent delivery system, an autosampler, and a diode-array detection (DAD). Zorbax Eclipse Plus C18 column (100 × 2.1 mm, 1.8 μm) was used for separation. The column temperature was set at 35 °C. The flow rate was 0.5 mL/min. The mobile phases were water with 0.1% formic acid (A) and acetonitrile with 0.1% formic acid (B), gradient elution was processed as table 1.

Table 1. Gradient elution program

时间(min)	A	B
0	91	9
9.5	91	9
25.5	83.5	16.5
33	77	23
44	46	54
45	0	100
50	0	100
50.1	91	9
60	91	9

Mass Spectrometry

Samples were analyzed by an Orbitrap Velos Pro mass spectrometer (Fig.1) using HESI- II probe in negative ion mode.

Ion source parameters: heat temperature, 300 °C; capillary temperature, 350 °C; spray voltage, 3.2 kV; sheath gas flow, 45 arb units; auxiliary gas flow, 5 arb units.

The resolution was set to 100,000 FWHM for full scan screening and 30,000 FWHM for Data dependent MSⁿ (DDA). Scan range was from 150 to 1400 amu, the isolation window for DDA acquisition was set as 2 amu, collision energy 40%.



Fig.1 Orbitrap Velos Pro mass spectrometer

Data Analysis

Data was processed by Thermo Scientific™ Compound Discoverer™ software for sulfur derivatives related peaks detection, paeoniflorin sulfonate was selected as reference molecule with the formula of C₂₃H₂₈O₁₃S (Fig.2), multiple nodes were constructed in a work flow to optimize the factors which may impact the precision of screening result.




Fig.2 Compound Discoverer software and interface of isotope pattern screening

Results

The separation of isotope pattern at ultra high resolution

Isotopic pattern filter was used mainly in study of Cl- and Br-containing compounds due to their distinct isotope patterns, but it has seldomly been used for screening of S-containing compounds, because the ^{34}S isotopic ion will be drastically affected by the contribution of $^{13}\text{C}_2$ and ^{18}O , which might leads to many false positives if the resolving power is not high enough. As the resolution capabilities have been increased along with the development of MS technique, mass spectrometers can now be able to separate the extremely close isobars and resolve the fine isotopic pattern of the small molecules. Take $\text{C}_{25}\text{H}_{23}\text{O}_{15}\text{S}$ for example, Figure 2 demonstrates the evolution of the A_2 as the resolution arrange from 30K to 100K. From a broad singlet to absolute two peaks, the improved resolution definitely enhances the recognition ability of the algorithm, which makes it possible to extract ^{34}S information directly.





Fig.3 Evolution of isotope pattern A_2 at different resolutions.

Sulfur derivatives search

There is no effective method for sulfur derivatives screening with conventional identification work flow, it will take lots of time to find them out from complex ms signal manually and may lead to confusing results: numerous formulas will be generated even the mass accuracy is better than 1ppm. Base on the ultra high resolution ms data provided by Orbitrap technology, Fine Isotope Pattern (FIP) filter can now be used to screen sulfur atom containing mass signal with specific tolerance (typically 5ppm mass tolerance for more accurate result and 60% intensity tolerance for more comprehensive result). significant sulfur-related peaks can be detected significantly after processed by FIP filter. The number of sulfur atoms can even be counted through the ratio of fine isotopes, then most of false positive formula could be excluded, leading to a more accurate result. Take 595.07672 m/z for example, ratio of ^{34}S is 4.6% and ratio of $^{13}\text{C}_2+^{18}\text{O}$ is 6.2% leading to the only element composition: $\text{C}_{25}\text{H}_{23}\text{O}_{15}\text{S}$.






Figure. 8. Comparison of FIP filter results at different mass tolerance

Conclusion

As the resolution capabilities have been increased along with the development of MS technique, mass spectrometers can now be able to resolve the fine isotope pattern, therefore, a novel fine isotopic pattern screening method was developed according to the separation of ^{34}S with $^{13}\text{C}_2$ and ^{18}O isotopic pattern which was successfully applied in study of sulfur-fumigated *Pueraria* species. The comparison result also shows that ultra high resolution, high mass accuracy and novel FIP algorithm are three major factors for this method.

Furthermore, the novel FIP method is also a potent tool for not only S-containing compounds but also other specific element such as ^{15}N or ^{18}O has a wide range of application:

- Discovery of specific element containing natural products from TCM
- Drug metabolites screening in complex biological matrix
- Definitive elemental composition determination in complex de novo applications

www.thermoscientific.com

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 10 3292 200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand